

Title:	How does the STX40 Shaker work			ID:
				0126
Date in:	Response:	Model:	Author:	
2006-08-31	2006-08-31	STX40	CMa	

**Q:**

Is there any information you could send me regarding the possible STX 44 shaker unit

**A:**

Max. amplitude is 6mm.

Up to 600 rpm the *maximum* amplitude is 3mm. If not otherwise specified we put in 1mm.

Max speed is 1'200 rpm. Speed can be controlled through the comport.

Shaker motion on access is controlled full automatically.