

Title:	Incubator with 260 MTP Capacity			ID:
				0146
Date in:	Response:	Model:	Author:	
2004-06-22	2004-06-22	STX500	C <i>Ma</i>	

Q:

We have another possible customer who needs to incubate 260 standard 96-well plates at a time.

Do you have a unit that large, or will he need several smaller incubators?

A:

We have the STX500 which you can have in three configurations:

STX500 (Tall stackers holding 53 MTP each)

STX502 (BiLevel 20x22 MTP standard stackers)

STX502 (BiLevel 10x22 MTP stackers on lower deck 10x28 MTP stackers on upper deck)

