

<b>Title:</b>	<b>Incubator Hotels</b>	<b>ID:</b>	
		0147	
<b>Date in:</b>	<b>Response:</b>	<b>Model:</b>	<b>Author:</b>
2003-01-08	2003-01-08	-	CMa

## Q:

I do actually have a few questions :

- 1) Do you have both incubators and hotels?
- 2) can the microplates be handled with lids?
- 3) If we are interested in your products is there a way to have a demo or a trial period?

## A:

- 1) Yes, we have both.
- 2) Yes, they can. The standard cassette (22 levels) is designed for MTP's with lid. If you want to be sure about the plates it is always a good idea to send sample plates for test.
- 3) Yes, we simply have to know the conditions, specifications the system has to meet