

<b>Title:</b>	<b>Power consumption of the STX40</b>			<b>ID:</b>
				0171
<b>Date in:</b>	<b>Response:</b>	<b>Model:</b>	<b>Author:</b>	
2004-04-27	2004-04-27	STX40	CMa	

**Q:**

How much is the maximum consumption of the incubator? (400W on the brochure, but 0.8KW on the tag on the back of the incubator.(See Tag.jpg))

**A:**

The power consumption varies a lot. On start-up the unit may draw as much as 800W. After a while it will average at 400W