


Title:	Shaker Data Information			ID:
				0181
Date in:	Response:	Model:	Author:	
2003-10-23	2003-10-23	-	CMa	

Q:

Do you offer a shaking incubator? Could you send me descriptions and costing information if you do offer this product.

A:

This is some information about the shaker. Basically it is an orbital shaker that can be remotely controlled. All control for the internal access – like positioning - is automatically done. So, the only thing you have to do once is start it and set it to the desired speed / frequency. The shaker is an option like the CO2. Frequency is <1..10 cycles per second. The default amplitude is 3mm pp and can be changed by replacement parts.



New StoreX Shaker

Agile Popular & Innovative

The advertisement features the LiCONiC INSTRUMENTS logo at the top right. It includes several images: a large view of the shaker's internal rollers, a close-up of a blue control panel with buttons, a red textured surface, and a view of the shaker's front door with a window.



New StoreX EasyClean

Simple Smart & Effective

The advertisement features the LiCONiC INSTRUMENTS logo at the top right. It includes several images: a view of the shaker's internal rollers, a view of the shaker's front door, and a woman in a white lab coat kneeling and cleaning the shaker with a spray bottle and a blue bucket.