


Title:	Shovel Motor pots for speed and power			ID:
				0190
Date in:	Response:	Model:	Author:	
2003-05-27	2003-05-27	-	CMA	

I need a picture of the shovel motor pots for speed and power, Some thing that points out which is which and how many turns is normal to change them. On the 3116 at 4C the shovel stalls sometimes.



Liconic Controller Card



Notes:

Some options may be missing depending on systems configuration
 CW (clockwise): Increasing
 CCW (counter clockwise): Decreasing
 Improper adjustments may lead to damage of equipment

STX_CtrlCard