

<b>Title:</b>	<b>Special Plates Height Cassettes</b>	<b>ID:</b>	
		0204	
<b>Date in:</b>	<b>Response:</b>	<b>Model:</b>	<b>Author:</b>
2004-06-21	2004-06-21	-	CMa

## Q:

I have just been given additional information about Novartis Boston, they will want to use 384 way Nunc plates, which are over 17mm high. Please let me know what you can do about this for us.

## A:

We have cassettes with different pitch. The next standard pitch is 28mm. We can also make special pitch cassettes. Note that your driver will have to send the configuration information for different pitches.