

Title:	STX500 Capacity	ID:	
		0216	
Date in:	Response:	Model:	Author:
2004-07-11	2004-07-11	STX500	CMa

Q:

A couple more StoreX 500 questions for you:

- 1) What is the actual capacity of the StoreX 500 with standard cassettes and how many plates are there per cassette?
- 2) What is the actual capacity of the StoreX 500 with the NUNC compatible cassettes used by TAP's Select-T and how many plates are there per cassette?

A:

- 1) 22 MTP height. You may use Standard stackers on the lower and CRS stackers on the upper deck. This will give you 500MTP capacity or take the tall stackers that give you 530 MTP capacity.
- 2) I strongly think there will be the same capacity. I'm expecting plate samples for confirmation.