

|                 |                            |               |                |
|-----------------|----------------------------|---------------|----------------|
| <b>Title:</b>   | <b>STX Family Overview</b> | <b>ID:</b>    |                |
|                 |                            | 0224          |                |
| <b>Date in:</b> | <b>Response:</b>           | <b>Model:</b> | <b>Author:</b> |
| 2003-02-19      | 2003-02-19                 | STX Series    | CMa            |

**Q:**

I'm interested in getting some basic information on your range of incubator products. I'm particularly interested in the robotic compatible versions with basic incubation and also with CO2. If you could forward me the features, specs and approximate US list pricing, I would greatly appreciate it.

**Q:**

Hope that the attached information will give you a first idea of the product range. The units are defined by

- Capacity (44..~1000 MTP)
- Climate Temperature (04..70°C)
- Climate Humidity (<5%RH or >90%RH)
- Configuration (StandAlone, BenchTop)

in addition there is a wide range of options (e.g. CO2, shaker ..) and accessories (e.g. stacker, transfer stations..).