

<b>Title:</b>	<b>Travel Length Slide Station</b>	<b>ID:</b>	
		0241	
<b>Date in:</b>	<b>Response:</b>	<b>Model:</b>	<b>Author:</b>
2005-05-13	2005-05-13	-	CMa

## Q:

Can you tell me what is the max length of the plate shuttle on a 9 tower plate hotel. I would like to put it at the back of an EVO & be able to access with RoMa.

## A:

The travel length you are using is 400mm. We can make almost any length. So far the longest was 1000mm