

<b>Title:</b>	<b>1401-2 vs 1701-2</b>			<b>ID:</b>
				316
<b>Date in:</b>	<b>Response:</b>	<b>Model:</b>	<b>Author:</b>	
2004-04-27	2004-04-27	STX40	CMa	

**Q:**

Is there such data that as of FW version on STR and STX series units we are now using ST 1401 and ST1402 verses ST 1701 ST 1702

**A:**

1701/2 is obsolete and is being replaced against 1401/2. 1401/2 on all V7x and STX44 from V6.9.07 and later. STR240 still have 1701/2