

Title:	LPR240-STR240 Worktable Hole Dimensions			ID:
				0377
Date in:	Response:	Model:	Author:	
2011-09-09	2004-04-27	LPR-240	CMA	

Q:

Can you tell me what the minimum size hole would be in a worktable to allow delivery of a MTP to a nest on the worktable? I looked at the housing for the LPR-240 and that hole is 260x270mm. However, I integrated an LPR-240 years ago with a hole that was 250x270mm so I know that the hole in the worktable can be smaller than the one in the housing. What I don't know is how much smaller it can be. Of course I am trying to preserve as many worktable grids as possible.

A:

If this is a new unit I recommend to put these "extras" into the order. On the opening we would have to make sure that the unit works with your dimensions. I remember that we added some room because the size of the opening was marginal.